

—A helping hand for your research

#### **Product Manual**

# Configuration-specific Monoclonal Antibody Based Arf6 Activation Assay Kit

1

Catalog Number: 82401

20 assays

# **Table of Content**

Product Description	3
Assay Principle	3
Kit Components	3
Storage	4
Materials Needed but Not Supplied	4
Reagent Preparation	4
Sample Preparation	5
Assay Procedure	6
Example of Results	8

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

NewEast Biosciences 24 Whitewoods Lane Malvern , PA 19355

### **Product Description**

ARF6 (ADP-ribosylation factor 6) is a member of the ARF super-family. ARF genes encode small GTPases that increase the ADP-ribosyltransferase activity of cholera toxin and are critical for vesicular trafficking as activators of phospholipase D. Arf6 regulates membrane trafficking and functions as a regulatory molecule of phagocytosis.

Currently there is no direct assay to measure the activation of Arf 6 GTPase.

NewEast Biosciences Arf6 Activation Assay Kit is based on the configuration-specific monoclonal antibody that specifically recognizes Arf6-GTP, but not Arf6-GDP, and an Arf6 specific rabbit polyclonal antibody. Given the high affinity of monoclonal antibodies to their antigens, the activation assay could be performed in a short time. This assay provides the reliable results with consistent reproducibility.

These anti-Arf 6-GTP monoclonal antibodies can also be used to monitor the activation of Arf 6 in cells and in tissues by immunohistochemistry.

NewEast Biosciences Arf6 Activation Assay Kit provides a simple and fast method to monitor the activation of Arf 6. Each kit provides sufficient quantities to perform 20 assays.

### Assay Principle

NewEast Biosciences Arf 6 Activation Assay Kit bases on the configuration-specific anti-Arf 6-GTP monoclonal antibody to measure the active Arf 6-GTP levels, either from cell extracts or from in vitro GTPγS loading Arf 6 activation assays. Briefly, anti-active Arf 6 mouse monoclonal antibody will be incubated with cell lysates containing Arf 6-GTP. The bound active Arf 6 will then be pulled down by protein A/G agarose. The precipitated active Arf 6 will be detected by immunoblot analysis using an anti-Arf 6 rabbit polyclonal antibody, respectively.

# Kit Components

- 1. <u>Anti-active Arf 6, Mouse Monoclonal Antibody (Catalog No. 26918):</u> One vial  $-22~\mu L$  (1 mg/ml) in PBS, pH 7.4, containing 50% glycerol and 0.05% sodium azide. This antibody specifically recognizes Arf 6-GTP from all vertebrates.
- 2. Protein A/G Agarose (Catalog No. 30301): One vial 400 μL of 50% slurry.

NewEast Biosciences 24 Whitewoods Lane Malvern, PA 19355

- 3. <u>5X Assay/Lysis Buffer (Catalog No. 30302):</u> One bottle 30 mL of 250 mM Tris-HCl, pH 8, 750 mM NaCl, 50 mM MgCl<sub>2</sub>, 5 mM EDTA, 5% Triton X-100.
- 4. <u>Anti- Arf 6, Rabbit Polyclonal Antibody (Catalog No. 21032):</u> One vial 100 μL (1 mg/ml) in PBS, pH 7.4, contained 50% glycerol.
- 5.  $\underline{100 \text{ X GTP}_{\gamma}\text{S}}$  (Catalog No. 30303): One vial -100  $\mu$ l at 10 mM, use 5  $\mu$ L of GTP $_{\gamma}\text{S}$  for GTP-labeling of 0.5 mL of cell lysate.
- 6. 100 X GDP (Catalog No. 30304): One vial -100 μl at 100 mM, use 5 μL of GDP for GDP-labeling of 0.5 mL of cell lysate.

# **Storage**

Store all kit components at 4°C until their expiration dates.

#### **Materials Needed but Not Supplied**

- 1. Stimulated and non-stimulated cell lysates
- 2. Protease inhibitors
- 3. 4 °C tube rocker or shaker
- 4. 0.5 M EDTA, pH8.0
- 5. 1 M MgCl<sub>2</sub>
- 6. 2X reducing SDS-PAGE sample buffer
- 7. Electrophoresis and immunoblotting systems
- 8. Immunoblotting wash buffer such as TBST (10 mM Tris-HCl, pH 7.4, 0.15 M NaCl, 0.05% Tween-20)
- 9. Immunoblotting blocking buffer (TBST containing 5% Non-fat Dry Milk or 3% BSA)
- 10. PVDF or nitrocellulose membrane
- 11. Secondary Antibody
- 12. ECL Detection Reagents

#### **Reagent Preparation**

• 1X Assay/Lysis Buffer: Mix the 5X Stock briefly and dilute to 1X in deionized water. Just prior to usage, add protease inhibitors such as 1 mM PMSF,  $10 \mu g/mL$  leupeptin, and  $10 \mu g/mL$  aprotinin.

NewEast Biosciences 24 Whitewoods Lane Malvern, PA 19355

# **Sample Preparation**

#### **Adherent Cells**

- 1. Culture cells (one 10-cm plate,  $\sim 10^7$  cells) to approximately 80-90% confluence. Stimulate cells with activator or inhibitor as desired.
- 2. Aspirate the culture media and wash twice with ice-cold PBS.
- 3. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cells (0.5 1 mL per 10 cm tissue culture plate).
- 4. Place the culture plates on ice for 10-20 minutes.
- 5. Detach the cells from the plates by scraping with a cell scraper.
- 6. Transfer the lysates to appropriate size tubes and place on ice.
- 7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3-4 times to shear the genomic DNA.
- 8. Clear the lysates by centrifugation for 10 minutes (12,000 x g at 4 °C).
- 9. Collect the supernatant and store samples (~1-2 mg of total proteins) on ice for immediate use, or snap freeze and store at 70 °C for future use.

#### **Suspension Cells**

- 1. Culture cells and stimulate with activator or inhibitor as desired.
- 2. Perform a cell count, and then pellet the cells by centrifugation.
- 3. Aspirate the culture media and wash twice with ice-cold PBS.
- 4. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cell pellet  $(0.5 1 \text{ mL per } 1 \text{ x } 10^7 \text{cells})$ .
- 5. Lyse the cells by repeated pipetting.
- 6. Transfer the lysates to appropriate size tubes and place on ice.
- 7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3-4 times to shear the genomic DNA.
- 8. Clear the lysates by centrifugation for 10 minutes (12,000 x g at 4 °C).
- 9. Collect the supernatant and store samples on ice for immediate use, or snap freeze and store at 70 °C for future use.

#### In vitro GTPyS/GDP Protein Loading for positive and negative controls

Note: In vivo stimulation of cells will activate approximately 10% of the available Arf 6, whereas in

NewEast BiosciencesFax: 610-945-200824 Whitewoods Lane5Phone: 610-945-2007Malvern , PA 19355Web: www.neweastbio.com

E-mail: sale@neweastbio.com

vitro GTPγS protein loading will activate nearly 90% of the Arf 6.

- 1, Aliquot 0.5 ml of each cell extract to two microfuge tubes (or use 1  $\mu$ g of purified Arf 6 protein).
- 2, To each tube, add 20 µl of 0.5 M EDTA (to 20 mM final concentration).
- 3, Add 5 μl of 100 X GTPγS (to 100 μM, final concentration) to one tube (positive control).
- 4, Add 5 µl of 100 X GDP (to 1 mM, final concentration) to the second tube (negative control).
- 5. Incubate the tubes at 30°C for 30 minutes with agitation.
- 6, Stop loading by placing the tubes on ice and adding  $32.5 \mu l$  of 1 M MgCl<sub>2</sub> (to 60 mM, final concentration).

#### **Assay Procedure**

#### I. Active Arf 6 Pull-Down Assay

- 1. Aliquot 0.5 1 mL of cell lysate to a microcentrifuge tube.
- 2. Adjust the volume of each sample to 1 mL with 1X Assay/Lysis Buffer.
- 3. Add 1 µl anti-active Arf 6 monoclonal antibody to the tube.
- 4. Thoroughly resuspend the protein A/G Agarose bead slurry by vortexing or titurating.
- 5. Quickly add 20 µL of resuspended bead slurry to each tube.
- 6. Incubate the tubes at 4 °C for 1 hour with gentle agitation.
- 7. Pellet the beads by centrifugation for 1 min at 5,000 x g.
- 8. Aspirate and discard the supernatant, making sure not to disturb/remove the bead pellet.
- 9. Wash the bead 3 times with 0.5 mL of 1X Assay/Lysis Buffer, centrifuging and aspirating each time.
- 10. After the last wash, pellet the beads and carefully remove all the supernatant.
- 11. Resuspend the bead pellet in 20 μL of 2X reducing SDS-PAGE sample buffer.
- 12. Boil each sample for 5 minutes.
- 13. Centrifuge each sample for 10 seconds at 5,000 x g.

#### II. Electrophoresis and Transfer

1. Load 15 μL/well of pull-down supernatant to a polyacrylamide gel (17%). Also, it's

NewEast Biosciences 24 Whitewoods Lane Malvern, PA 19355


- recommended to include a pre-stained MW standard (as an indicator of a successful transfer in step 3).
- 2. Perform SDS-PAGE following the manufacturer's instructions.
- 3. Transfer the gel proteins to a PVDF or nitrocellulose membrane following the manufacturer's instructions

#### III. Immunoblotting and Detection (all steps are at room temperature, with agitation)

- 1. Following the electroblotting step, immerse the PVDF membrane in 100% Methanol for 15 seconds, and then allow it to dry at room temperature for 5 minutes.
  - *Note: If Nitrocellulose is used instead of PVDF, this step should be skipped.*
- 2. Block the membrane with 5% non-fat dry milk or 3% BSA in TBST for 1 hr at room temperature with constant agitation.
  - Incubate the membrane with anti- Arf 6 rabbit polyclonal antibody, freshly diluted 1:50~1000 (depending on the amount of Arf 6 proteins in your samples) in 5% non-fat dry milk or 3% BSA/TBST, for 1-2 hr at room temperature with constant agitation or at 4°C overnight.
- 3. Wash the blotted membrane three times with TBST, 5 minutes each time.
- 4. Incubate the membrane with a secondary antibody (e.g. Goat Anti-Rabbit IgG, HRP-conjugate), freshly diluted 1:1000 in 5% non-fat dry milk or 3% BSA/TBST, for 1 hr at room temperature with constant agitation.
- 5. Wash the blotted membrane three times with TBST, 5 minutes each time.
- 6. Use the detection method of your choice such as ECL.

# **Example of Results**

The following figure demonstrates typical results seen with NewEast Biosciences Arf 6 Activation Assay Kit. One should use the data below for reference only.



IP: anti-active Arf6 mAb

IB: anti-Arf6 pAb

**Arf6 activation assay.** Purified Arf6 proteins were immunoprecipitated after treated with GDP (lane 1) or GTPγS (lane 2). Immunoprecipitation was done with the anti-active Arf6 monoclonal antibody (Cat. No. 26918). Immunoblot was with an anti-Arf6 rabbit polyclonal antibody(Cat. No. 21032).