



—A helping hand for your research

---

Configuration-specific Monoclonal Antibody Based

## **Gα13 Activation Assay Kit**

(30 Assays)

**Cat. # 80401**

---

Address: NewEast Biosciences  
24 Whitewoods Lane  
Malvern, PA 19355

Orders: 610-945-2007  
sale@neweastbio.com  
Support: 610-945-2007  
info@neweastbio.com  
Web: www.neweastbio.com

# Gα13 Activation Assay Kit Protocol

Cat. # 80401

**FOR RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC PROCEDURES**

## Product Description

A structurally diverse repertoire of ligands, from photons to large peptides, activates GPCRs to elicit their physiological functions. Ligand-bound GPCRs, in turn, function as guanine nucleotide exchange factors catalyzing the exchange of GDP bound on the Gα subunit with GTP in the presence of Gβγ, causing the dissociation of the Gα subunit from the Gβγ dimer to form two functional units (Gα and Gβγ). Both Gα and Gβγ subunits signal to various cellular signaling pathways. Based on the sequence and functional homologies, G proteins are grouped into four families: Gs, Gi, Gq, and G12. As increasing numbers of effectors and interacting proteins for these G proteins have been identified, the physiological processes in which G proteins participate are multiplying.

Among the four subfamilies of G proteins, the function of G12/13 subfamily is less well understood. In this family, there are two members, G12 and G13, that are expressed ubiquitously. Gα12 knockout mice appeared normal. Gα13 knockout mice displayed embryonic lethality (~E9.5). The Gα13 mouse embryos had defective vascular systems. Gα13 is also essential for receptor tyrosine kinase-induced migration of fibroblast and endothelial cells.

NewEast Biosciences Gα13 Activation Assay Kit provides a simple and fast tool to monitor the activation of Gα13. Each kit provides sufficient quantities to perform 30 assays.

## Assay Principle

NewEast Biosciences Gα13 Activation Assay Kit bases on the configuration-specific-anti-Gα13-GTP monoclonal antibody to measure the active Gα13-GTP levels, either from cell extracts or from in vitro GTPγS loading Gα13 activation assays. Briefly, antiactive Gα13 mouse monoclonal antibody will be incubated with cell lysates containing Gα13-GTP. The bound active Gα13 will then be pulled down by protein A/G agarose. The precipitated active Gα13 will be detected by immunoblot analysis using anti-Gα13 rabbit polyclonal antibody.

## Kit Contents


This kit contains enough reagents for approximately 30–35 pull-down assays.

Reagent	Cat. #	Quantity	Storage
Anti – active Gα13 Mouse Monoclonal Antibody	Cat. # 26902	1 X 35μl	-20°C
Protein A/G Agarose	Cat. # 30301	1X600μl	4°C
5X Assay/Lysis Buffer	Cat. # 30303	1X30mL	4°C
Anti-Gα13 Rabbit polyclonal Antibody	Cat. # 21005	1X50μl	-20°C
100x GTPγS	Cat. # 30302	1X50μl	--80°C
100x GDP	Cat. # 30304	1X50μl	-80°C
HRP – Goat Anti-Rabbit IgG	Cat. # 29002	1X50μl	-20°C

Note: For GDP and GTPγS, aliquot into 10x5ul volumes, then store at -80 degrees.

## Example of Results

The following figure demonstrates typical results seen with NewEast Biosciences Gα13 Activation Assay Kit. One should use the data below for reference only.



**Gα13 activation assay.** MEF cells were treated with (lane 2) or without (lane 1) LPA. Cell lysates were incubated with an antiactive Gα13 monoclonal antibody (Cat. # 26902) (top panel). The precipitated active Gα13 was immunoblotted with an anti-Gα13 rabbit polyclonal antibody (Cat # 21005). The bottom panel shows the Western blot with antiGα13 of the cell lysates used (5% of that used in the top panel).

# Gα13 Activation Assay Kit Protocol

Cat. # 80401

## Materials Needed but Not Supplied

Stimulated and non-stimulated cell lysates  
Protease inhibitors  
4°C tube rocker or shaker  
0.5 M EDTA, pH8.0  
1 M MgCl<sub>2</sub>  
2X reducing SDS-PAGE sample buffer  
Electrophoresis and immunoblotting systems  
Immunoblotting wash buffer such as TBST  
(10 mM Tris-HCl, pH 7.4, 0.15 M NaCl, 0.05% Tween-20)  
Immunoblotting blocking buffer  
(TBST containing 5% Non-fat Dry Milk or 3% BSA)  
ECL Detection Reagents

## A Reagent Preparation

1X Assay/Lysis Buffer: Mix the 5X Stock briefly and dilute to 1X in deionized water. Just prior to usage, add protease inhibitors such as 1 mM PMSF, 10 µg/mL leupeptin, and 10 µg/mL aprotinin

## B Sample Preparation

### Adherent Cells

1. Culture cells (one 10-cm plate, ~ 10<sup>7</sup> cells) to approximately 80–90% confluence. Stimulate cells with activator or inhibitor as desired.
2. Aspirate the culture media and wash twice with ice-cold PBS.
3. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cells (0.5– 1 mL per 10 cm tissue culture plate).
4. Place the culture plates on ice for 10–20 minutes.
5. Detach the cells from the plates by scraping with a cell scraper.
6. Transfer the lysates to appropriate size tubes and place on ice.
7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3–4 times to shear the genomic DNA.
8. Clear the lysates by centrifugation for 10 minutes (12,000 x g at 4°C).
9. Collect the supernatant and store samples (~1–2 mg of total proteins) on ice for immediate use, or snap freeze and store at – 70°C for future use.

### Suspension Cells

1. Culture cells and stimulate with activator or inhibitor as desired.
2. Perform a cell count, and then pellet the cells by centrifugation.
3. Aspirate the culture media and wash twice with ice-cold PBS.
4. Completely remove the final PBS wash and add ice-cold 1X Assay/Lysis Buffer to the cell pellet (0.5 – 1 mL per 1 x 10<sup>7</sup> cells).
5. Lyse the cells by repeated pipetting.
6. Transfer the lysates to appropriate size tubes and place on ice.
7. If nuclear lysis occurs, the cell lysates may become very viscous and difficult to pipette. If this occurs, lysates can be passed through a 27½-gauge syringe needle 3–4 times to shear the
8. Clear the lysates by centrifugation for 10 minutes (12,000 xg at 4°C).
9. Collect the supernatant and store samples on ice for immediate use, or snap freeze and store at – 70°C for future use.

## C In vitro GTP γ S/GDP Protein Loading for positive and negative controls

Note: In vivo stimulation of cells will activate approximately

10% of the available Gα13, whereas in vitro GTP γ S protein loading will activate nearly 90% of Gα13.

1. Aliquot 0.5 ml of each cell extract to two microfuge tubes (or use 1 µg of purified Gα13 protein).
2. To each tube, add 20 µl of 0.5 M EDTA (to 20 mM final concentration).
3. Add 5 µl of 100 X GTP γ S (to 100 µM, final concentration) to one tube (positive control).
4. Add 5 µl of 100 X GDP (to 1 mM, final concentration) to the second tube (negative control).
5. Incubate the tubes at 30°C for 30 minutes with agitation.
6. Stop loading by placing the tubes on ice and adding 32.5 µl of 1 M MgCl<sub>2</sub> (to 60 mM, final concentration).

## D Affinity Precipitation of Activated G protein

1. Aliquot 0.5 – 1 mL of cell lysate (~1 mg of total cellular protein) to a microcentrifuge tube.
2. Adjust the volume of each sample to 1 mL with 1X Assay/Lysis Buffer.
3. Add 1 µl anti-active Gα13 monoclonal antibody to the tube.
4. Thoroughly resuspend the protein A/G Agarose bead slurry by vortexing or titrating.
5. Quickly add 20 µL of resuspended bead slurry to each tube.
6. Incubate the tubes at 4°C for 1 hour with gentle agitation.
7. Pellet the beads by centrifugation for 1 min at 5,000 x g.
8. Aspirate and discard the supernatant, making sure not to disturb/remove the bead pellet.
9. Wash the bead 3 times with 0.5 mL of 1X Assay/Lysis Buffer, centrifuging and aspirating each time.
10. After the last wash, pellet the beads and carefully remove all the supernatant.
11. Resuspend the bead pellet in 20 µL of 2X reducing SDS-PAGE sample buffer.
12. Boil each sample for 5 minutes.
13. Centrifuge each sample for 10 seconds at 5,000 x g

## E Western blot analysis

1. Load 15 µL/well of pull-down supernatant to a polyacrylamide gel (17%). Also, it's recommended to include a pre-stained MW standard (as an indicator of a successful transfer in step 3).
2. Perform SDS-PAGE following the manufacturer's instructions.
3. Transfer the gel proteins to a PVDF or nitrocellulose membrane following the manufacturer's instructions
4. Following the electroblotting step, immerse the PVDF membrane in 100% Methanol for 15 seconds, and then allow it to dry at room temperature for 5 minutes.  
*Note: If Nitrocellulose is used instead of PVDF, this step should be skipped.*
5. Block the membrane with 5% non-fat dry milk or 3% BSA in TBST for 1 hr at room temperature with constant agitation. Incubate the membrane with anti-Gα13 polyclonal antibody, freshly diluted 1:50–500 (depending on the amount of Gα13 proteins in your samples) in 5% non-fat dry milk or 3% BSA/TBST, for 1–2 hr at room temperature with constant agitation or at 4°C overnight.
6. Wash the blotted membrane three times with TBST, 5 minutes each time.
7. Incubate the membrane with a secondary antibody (e.g. Goat Anti-Rabbit IgG, HRP-conjugate), freshly diluted 1:1000 in 5% non-fat dry milk or 3% BSA/TBST, for 1 hr at room temperature with constant agitation.
8. Wash the blotted membrane three times with TBST, 5 minutes each time.
9. Use the detection method of your choice such as ECL