

## **Product Description**

A helping hand for your research

## Ral

Catalog Number: 21150 Gene Symbol: RAL

**Description:** Anti-Ral Rabbit Polyclonal Antibody

Background: The Ras-like small G proteins, RalA/B, are important components of Ras signaling pathways, implicated in the initiation and maintenance of tumorigenic transformation, vesicle as well as transport, apoptosis, transcription, cell migration,

proliferation.

Immunogen: Two synthetic peptides from the internal regions of RalA and RalB proteins. Tested applications: ELISA, WB, IHC, IP

## **Recommended dilutions:**

ELISA 1:100-1:5000 WB 1:100-1:2000 IHC 1:100-1:1000 Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal Purity: Purified from serum

Format: Liquid Storage buffer: Preservative: no

Constituents: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH

7.4, 150 mM NaCl, 50% glycerol


Species Reactivity: Ral antibody recognizes

RalA and RalB from vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze

/ thaw cycles

## Western blot:



EC-1 cell lysate was used for Western blot with anti-RalA polyclonal antibody (cat # 21150).

Web: www.neweastbio.com