


## RalA

Catalog Number: 21034 Gene Symbol: RALA Description: Anti-RalA Rabbit Polyclonal Antibody

Background: The Ras-like small G proteins, RalA/B, are important components of Ras signaling pathways, implicated in the initiation and maintenance of tumorigenic transformation, well vesicle as transport, apoptosis, as transcription. cell migration, and cell proliferation. Immunogen: A synthetic peptide from the internal region of RalA protein Tested applications: ELISA, WB, IHC **Recommended dilutions:** ELISA 1:100-1:5000 WB 1:100-1:2000 IHC 1:100-1:1000 Concentration: 0.2 mg/ml Host: Rabbit Clonality: Polyclonal Purity: Purified from serum Format: Liquid Storage buffer: Preservative: no Constituents: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol Species Reactivity: RalA antibody recognizes RalA from vertebrates. Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



Western blot analysis of RalA proteins. HEK-293T cell lysate was used for Western blot with anti-RalA polyclonal antibody (cat # 21034).

For research use only. Not for diagnostic or therapeutic applications.