



Product Description

— A helping hand for your research

Warning:

This product is for research only, not for human or veterinary diagnostic or therapeutic use !

Product Name : Rabbit anti-FITC polyclonal antibody

Catalog NO. : 20010

Immunogen : KLH-FITC

Isotype : IgG

Size : 1 mg

Purification Method : Affinity-purified

Purity : >98% by SDS-PAGE

Applications : ELISA(1.5-2.5 μ g/mL Recommended), IHC(0.5-20 μ g/mL Recommended)

Buffer : 0.01 M PBS,50% glycerol with 0.05% Proclin 300

Storage : Store at -20°C for three years and avoid freeze/thaw cycles

24 Whitewoods Lane, Malvern, PA 19355 Tel: (610) 945-2007 Fax: (610) 945-2008

E-mail:sale@neweastbio.com Web:www.neweastbio.com