

Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

Rab7

Catalog Number: 26055 Gene Symbol: RAB7 Description: Anti-Rab7 Mouse Monoclonal Antibody

Background: Rab7 belongs to the Rab family of Western blot: small GTPases. Rab7 is located in the late endosome, lysosome, cytoplasmic vesicle, - 72 phagosome and melanosome. It is involved in the late endocytic transport and contributes to the - 55 maturation of phagosomes. Immunogen: purified His-Rab7 protein Tested applications: ELISA, WB, IHC, IP - 36 **Recommended dilutions:** 28 ELISA 1:100-1:5000 WB 1:100-1:2000 1:100-1000 IHC Concentration: 0.8 mg/ml - 17 Host: Mouse Clonality: Monoclonal Isotype: IgG Western blot analysis of recombinant Rab7 Purity: Purified from ascites proteins. Purified His-tagged Rab7 protein (Cat. Format: Liquid #10167) was blotted with anti-Rab7 monoclonal Storage buffer: antibody (Cat # 26055). Preservative: no Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol Species Reactivity: Rab7 antibody recognizes Rab7 proteins of vertebrates. Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles.