


Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

Rab11

Catalog Number: 26030 Gene Symbol: RAB11 Description: Anti-Rab11 Mouse Monoclonal Antibody

Background: Small GTPases of the Rab-family are crucial elements of the machinery that controls membrane traffic. Rab11 has been shown to control traffic through the recycling endosome. Immunogen: Recombinant full-length protein of Rab11A. Tested applications: ELISA, WB, IHC, IP **Recommended dilutions:** ELISA 1:100-1:5000 WB 1:100-1:1000 IHC 1:50-1:100 **Concentration:** 1 mg/ml Host: mouse **Clonality:** Monoclonal Isotype: IgG Purity: Purified from ascites Format: Liquid Storage buffer: Preservative: no Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol Species Reactivity: Rab11 antibody recognizes Rab11 proteins of vertebrates. Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



WB: anti-Rab11 mAb

Western blot analysis of recombinant Rab11A proteins (Cat# 26030).