

Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

RAB35


Catalog Number: 26033 Gene Symbol: RAB35 Description: Anti-Rab35 Mouse Monoclonal Antibody

Background: The Rab family of GTPases regulates eukaryotic vesicular membrane traffic. Rab35 is with both plasma membrane and endosomal localization, and has been implicated in diverse processes that include T-cell receptor recycling, oocyte yolk protein recycling and cytokinesis. Rab35 regulates neurite outgrowth in neuronal-like cells, and can induce protrusions.

Immunogen: Recombinant full length protein of Rab35.

Tested applications: ELISA, WB, IHC, IP **Recommended dilutions:** ELISA 1:100-1:5000 WB 1:100-1:2000 IHC 1:100-1:1000 Concentration: 1 mg/ml Host: mouse **Clonality:** Monoclonal Isotype: IgG Purity: Purified from ascites Format: Liquid Storage buffer: Preservative: no Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol Species Reactivity: Rab35 antibody recognizes Rab35 of vertebrates. Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



WB:anti-Rab35 mAb

Western blot analysis of recombinant GST-Rab35 proteins.