Wuhan NewEast Biosciences Co.Ltd



Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

Ras

Catalog Number: 26014 Gene Symbol: RAS

Description: Anti-Ras Mouse Monoclonal Antibody

Background: Small GTPases are a super-family of cellular signaling regulators. Ras belongs to the Ras sub-family of GTPases that regulate cell proliferation, cell motility, and gene transcription. GTP binding increases the activity of Ras, and the hydrolysis of GTP to GDP renders it inactive. GTP hydrolysis is aided by GTPase activating proteins (GAPs), while exchange of GDP for GTP is facilitated by guanine nucleotide exchange factors (GEFs).

Immunogen: Recombinant full length protein of

Ras

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA 1:100-1:5000 WB 1:100-1:2000 Concentration: 1 mg/ml

Host: mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid **Storage buffer:** Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH

7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Ras antibody recognizes H-Ras, N-Ras and K-Ras from vertebrates. **Storage Conditions:** Store at -20°C. Avoid

freeze / thaw cycles

Western blot:



WB: anti-Ras mAb

Western blot analysis of Ras proteins. MEF whole-cell lysate was used for western blot with the monoclonal anti-Ras antibody (cat # 26014).