Wuhan NewEast Biosciences Co.Ltd



Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

RhoA

Catalog Number: 26007 Gene Symbol: RHOA

Description: Anti-RhoA Mouse Monoclonal Antibody

Background: Small GTPases are a super-family of cellular signaling regulators. RhoA belongs to the Rho sub-family of GTPases. Rho proteins play critical roles in many actin cytoskeleton-dependent processes including platelet aggregation, cell motility, contraction, and cytokinesis. It regulates the formation of stress fibers and focal adhesions in fibroblasts and Ca²⁺ sensitivity of smooth muscle contraction.

Immunogen: Recombinant full length protein of

RhoA.

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA 1:1000-1:5000 WB 1:1000-1:2000

Concentration: 1 mg/ml

Host: mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid **Storage buffer:** Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH

7.4, 150 mM NaCl, 50% glycerol


Species Reactivity: RhoA antibody recognizes

RhoA from vertebrates.

Storage Conditions: Store at -20°C. Avoid

freeze / thaw cycles

Western blot:



WB: anti-RhoA mAb

Western blot analysis of RhoA proteins. MEF whole-cell lysate was used for western blot with the anti-RhoA mouse monoclonal antibody (Cat. # 26007).