


Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

## **Apolipoprotein A1**

Catalog Number: 26025 Gene Symbol: APOA1 Description: Anti-Apolipoprotein A1 Mouse Monoclonal Antibody

Background: Apolipoprotein A1 is the major protein component of high-density lipoprotein (HDL) ("good cholesterol") in plasma. Apolipoprotein A1 helps to clear cholesterol from arteries. Immunogen: Recombinant full length protein of human Apolipoprotein A1 Tested applications: ELISA, WB, IHC, IP **Recommended dilutions:** ELISA 1:100-1:5000 WB 1:100-1:2000 IHC 1:50-1:100 Concentration: 1 mg/ml Host: mouse Clonality: Monoclonal Isotype: IgG Purity: Purified from ascites Format: Liquid Storage buffer: Preservative: no Constituents: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol Species Reactivity: Apolipoprotein A1 antibody recognizes Apolipoprotein A1 of primates. Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

## Western blot:



WB: anti-Apo A1 mAb

Western blot analysis of Apolipoprotein A1 proteins.