



Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

SMAD4

Catalog Number: 26453

Gene Symbol: SMAD4, DPC4, JIP, MADH4, MYHRS **Description:** Anti-SMAD4 Mouse Monoclonal Antibody

Background: SMAD4 family member 4, also known as SMAD4, is a protein that in human is encoded by the SMAD4 gene. It is often found mutated in many cancers. It acts as a tumor suppressor that functions in the regulation of the TGF-β signal transduction pathway, which negatively regulates growth of epithelial cells and the extracellular matrix (ECM). SMAD4 alterations have been found in multiploid colorectal cancer and pancreatic carcinoma. It is found inactivated in at least 50% of pancreatic cancers. It is also found mutated in the autosomal dominant disease juvenile polyposis syndrome (JPS).

Immunogen: Purified SMAD4 protein, human

origin.

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA: 1:1000-1:2000 WB: 1:500-1:1000 IHC: 1:50-1:100 Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal **Purity:** Purified from ascites

Format: Liquid Storage buffer: Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH

7.4, 150 mM NaCl, 50% glycerol


Species Reactivity: Recognizes SMAD4 of

vertebrates.

Storage Conditions: Store at -20°C. Avoid

freeze / thaw cycles.

Western blot:




Western blot analysis of recombinant SMAD4 protein. Purified His-tagged SMAD4 protein (1-552aa, 65kd) was blotted with anti-SMAD4 monoclonal antibody (Cat. #26453).





Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

Immunofluorescence:



Immunofluorescence of cells expressing SMAD4 proteins with anti-SMAD4 antibody.

HEK293T cells were transfected with pCDNA3-GFP-SMAD4 (WT) plasmid, then fixed and stained with anti-SMAD4 monoclonal antibody (Cat. #26453).