

Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

$G\pmb{\alpha}_{\!s}$

Catalog Number: 21007 **Gene Symbol:** Gnas **Description:** Anti-Gα_s Rabbit Polyclonal Antibody

Background: Heterotrimeric G proteins are critical cellular signal transducers. $G\alpha_s$ represents		Western blot:
one sub-family of G proteins that could mediate		
the stimulation of adenylyl cyclases.		
Immunogen: Recombinant full length protein of		MEF cells
Gα _s		kD
Tested applications: ELISA, WB, IHC		72 –
Recommended dilutions:		55 –
ELISA	1:100-1:5000	55 -
WB	1:100-1:2000	43 – — – G α _s
Concentration: 1 mg/ml		33 -
Host: Rabbit		33 -
Clonality: Polyclonal		26 –
Purity: Purified from serum		WB: anti-Gα _S pAb
Format: Liquid		
Storage buffer:		Western blot analysis of $G\alpha_s$ proteins. MEF
Preservative: no		whole-cell lysate was used for western blots.
Constituents: PBS (without Mg ²⁺ and Ca ²⁺), pH		
7.4, 150 mM NaCl, 50% glycerol		
Species Reactivity: $G\alpha_s$ antibody recognizes		
$G\alpha_s$ of vertebrates		

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles