

Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

Sar1a

Catalog Number: 21090 Gene Symbol: SAR1A Description: Anti-Sar1a Rabbit Polyclonal Antibody

Background:Small GTPases are a super-family
of cellular signaling regulators. Sar1 is a member
of these small GTPases. Sar1 functions as a
molecular switch to control protein—protein and
protein—lipid interactions that direct vesicle
budding from the ER.Immunogen:A synthetic peptide from the
internal region of human Sar1a proteinTested applications:ELISA, WB, IHCRecommended dilutions:ELISAELISA1:1000-1:5000

 WB
 1:100-1:2000

 IHC
 1:100-1:1000

 Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal

Purity: Purified from antisera

Format: Liquid

Storage buffer:

Preservative: no


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH

7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Sar1a antibody recognizes Sar1a proteins of vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles





Western blot analysis of Sar1a proteins. MEF whole-cell lysates were used for Western blot with anti-Sar1a rabbit polyclonal antibody (Cat# 21090).