Wuhan NewEast Biosciences Co.Ltd



Email:neweastbio@126.com Web:www.neweastbio.com.cn Tel:027-87561033

RalB

Catalog Number: 21033 Gene Symbol: RALB

Description: Anti-RalB Rabbit Polyclonal Antibody

Background: The Ras-like small G proteins, RalA/B, are important components of Ras signaling pathways, implicated in the initiation and maintenance of tumorigenic transformation, as well as vesicle transport, apoptosis, transcription, cell migration, and cell proliferation.

Immunogen: A synthetic peptide from the

internal region of RalB protein.

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA 1:100-1:5000 WB 1:100-1:2000 IHC 1:100-1:1000 Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal

Purity: Purified from serum

Format: Liquid Storage buffer: Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH

7.4, 150 mM NaCl, 50% glycerol


Species Reactivity: RalB antibody recognizes

RalB from vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze

/ thaw cycles

Western blot:



WB: anti-RalB pAb

Western blot analysis of RalB proteins. EC-1 cell lysate was used for Western blot with anti-RalB

polyclonal antibody (cat # 21033).